


Helios Occlusion Balloon Catheter

The Helios occlusion balloon catheter is specifically designed to work with the ImageWire, offering a safe method to clear blood from the viewing field. The Helios nominal inflation pressure is only 0.3 to 0.5 ATM, with no vessel distortion or vessel wall damage.



Specifications

Rated Burst Pressure:

1ATM minimum rated burst pressure

Balloon Length:

2.0 mm

Balloon Diameter:

4.0 mm (nominal)

Crossing Profile:

≤ 4.3 F (wrapped balloon)